



Flexible positioning.
Please note: The chair has a 25 stone weight capacity but most carers will struggle to manoeuvre clients over 15 stone in weight.



4 large swivel castors, 2 at rear are lockable.

Model shown in Zest Lime with black VP cushioning and visco foam quilt.



Optional profile head cushion



Full recline position



High leg lift position



High rise position

Porter optional extras



Legrest extension



Removable footplate length adjustable



Removable footplate Length & angle adjustable




Adjustable sliding, lockable thoracic supports (pair)



Soft fibre filled lateral supports (pair)



Careform backrest (Tredegar style)



Wedge back



Butterfly back



'Aircell' adjustable pressure waterfall back cushions




Primagel Gel pressure relief (medium / high risk)



Freedom pack (Accupack) battery back up 40 cycles



Head Cushion



Profile head cushion



Lap strap

Sizes

Porter chairs can be made in many sizes. Typical seat heights are a minimum 18".

A standard electric Single tilt-in-space Porter is:-

Seat height:- 18"
Seat depth:- 20"
Seat width:- 20"
Overall chair width:- 28"

Your local primacare distributor is:-

tel: 01685 845900 fax: 01685 844538
www.primacare.co.uk info@primacare.co.uk



The Porter Range



BRITISH DESIGN - BRITISH MANUFACTURE

The Porter Range



Foot tray extended. Seat can be extended forwards. TV position. Relaxed sleeping position.

Push back manual

1. Easy to use manual recliner Porter chair
2. Extra long extending legrest.
3. Removable seat module including Visco seat foam topper as standard.
4. Push handle for Carer to push chair or pull back to recline chair.
5. Sliding foot tray
6. 4 swivel castors, 2 at rear braked.



Adjusted to a fully upright position. Reclined. Tilted back. Extra back angle recline.

Dual tilt-in-space

Electric & Manual versions

1. Smooth electric tilt-in-space recline action, or manual legrest operated by finger pull handle and handle with trigger for tilt-in-space recline.
2. The dual tilts in space in an identical manner to the single motor chairs. The extra motor allows adjustable back angle positioning.
3. Sliding foot tray.
4. Removable cushioning for seat depth alteration. Continuous leg support quilt with Visco foam as standard.
5. Four swivel castors, 2 at rear braked.



Foot tray extended. Lateral support cushions (optional extra)



Normal recline. High leg recline.



Fully removable seat and quilt. Different seat cushion modules available

Single tilt-in-space

Electric & Manual versions

1. Smooth electric tilt-in-space recline action, or manual legrest operated by finger pull handle and handle with trigger for tilt-in-space recline.
2. Sliding foot tray.
3. Removable cushioning for seat depth alteration. Continuous leg support quilt with Visco foam as standard.
4. Four swivel castors, 2 at rear braked.



Bariatric BLTR 40 stone

Primacare can make porter chairs to whatever size and action type are required.

The chair shown here is a 40 stone heavy leg lift (20 stone) porter chair with a seat width of 30". The chair and wheels are designed to support occupants but this would be in a static situation.

The mobile aspect of this chair is to move the chair without the service user in the chair. 2 nurses would not be able to push an occupant with a weight of 40 stone successfully. An optional electrically powered 5th wheel can be added to assist when portering heavy occupants.



Chair tilted with high extending leg rest. Strong leg rest (20 stone). Power assisted portering.

Infection control

Some care homes and hospitals are now putting more emphasis on features of the chair that aid cleaning and do not harbour traps for contamination.

All Push back, Single and Dual chairs can be manufactured to the infection control specification: Fully PVC covered areas for cleaning using a mild bleach solution, and no velcro or zips that can trap dirt. All the cushioning is attached using magnets which enable all the seat, quilt and back cushions to be removed quickly for cleaning.

Cushions quickly removed for easy cleaning



Hoist friendly

This version has a raised seat height to enable access under the chair if required for hoist supports. The clearance is 5" (125mm).

All Push back, Single and Dual chairs can be manufactured to the Hoist friendly specification

