

Affinity

Bariatric Porter

The Affinity Bariatric Porter has a weight capacity of 40 or 50 stone and uses the reliable BLTR actions.


This gives the user separate control of the Back angle, Leg rest, Tilt and Rise functions of the porter chair.

With chairs of this weight capacity the portering aspect is to allow the chair to be moved around without the user in the chair. It would be very difficult for 2 carers to move the chair with a 50 stone occupant



primacare | always there

Affinity *Bariatric Porter* | specification.



		Standard
A	seat height	16-22"
B	seat width	24-34"
C	seat depth	18-26"
D	back height	27-33"
E	back angle	40° plus 30° tilt
F	overall width	Seat plus 9"
G	overall depth	38"
H	overall height	Seat height plus back height minus 3"

Standard features

- BLTR action with separate control of Back, Leg, Tilt and Rise functions
- Ultima vinyl Crib 5 with waterproof vapour permeable 4 way stretch fabric on all contact surfaces
- Removable seat quilt system giving access to the seat foam block.
- 4" Dia swivel castors with braked rear castors
- 25 stone extra heavy duty leg rest
- Flat vertical lift if required

10

Affinity *Bariatric Porter* | features and options.



The Bariatric Porter can be reclined and tilted to create an almost bed like flat position



This is a heavy chair and one career will only move the chair small distances with the user seated.



The Porter chair has a good high get out position



The chair can be raised 4" vertically